BeautiLink SA

- en** Instructions for use
- es** Instrucciones de uso
- fr** Mode d'emploi
- pt** Instruções de uso

SHOFU INC.
11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan
www.shofu.com

SHOFU DENTAL CORPORATION
1225 Stone Drive, San Marcos, California 92078-4059, USA
SHOFU DENTAL ASIA-PACIFIC PTE.LTD.
1 Science Park Road, #01-01/02 The Capricorn, Singapore Science Park II, Singapore 117528

IFU-000003
01.0.2023-10
72940

en READ CAREFULLY BEFORE USE

Self-adhesive resin cement

PRODUCT DESCRIPTION

BeautiLink SA is a dual-curing, self-adhesive resin cement for alumina, zirconia, silicate glass (lithium disilicate etc.), metal, and composite resin restorations. It has a choice of automix delivery or handmix. To enhance adhesion of indirect restorations with low retentive surfaces (crowns, table top onlays etc.) to the tooth substance, BeautiBond Xtreme can be used in conjunction with this material. BeautiLink SA is available in the three shades, Clear, Ivory and Opaque. Mixing ratio of the two pastes must be 1vol/1vol.

INTENDED PURPOSE

(Intended purpose)

BeautiLink SA is intended for cementation of indirect restorations.

(Intended patient population)

Patients with permanent teeth

(Clinical benefit)

To set the indirect restorations in place

INDICATIONS FOR USE

- Cementation of all-ceramic, composite and metal restorations
- Post cementation

PRECAUTIONS AND WARNINGS

1. Do not use this product if a patient or an operator is known to be allergic to this product or any of the material's ingredients.
2. If any inflammation or other allergic reactions occur on either the patient or operator, immediately discontinue use and seek medical advice.
3. Use medical (dental) gloves and protective glasses to avoid sensitization to this product. In case of accidental contact with oral soft tissue or skin, immediately blot with alcohol moistened cotton ball, and rinse with plenty of water. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
4. In case any serious health hazards occur on a patient or an operator by use of this product, please report to local competent authority, our local distributor and sales company. Contact information on our sales companies is available on www.shofu.com.

DIRECTIONS FOR USE


Preparation

Allow the materials to reach room temperature before use when stored in the refrigerator.

- In cases ambient temperature is higher than 23 °C/73 °F, working time will be shortened, while ambient temperature is lower than 23 °C/73 °F, working time will be extended.

Preparation of the automix syringe

1. Turn the cap a 1/4 turn counterclockwise to align the projections on the caps with the grooves of the syringe. Remove the cap by holding the foot end of the cap and pulling downward. (Fig. A)
2. Depress the plunger to dispense the materials to confirm that the two pastes are flowing evenly at the same level. (Fig. B)
 - Check this before every use.
3. For cementation of restorations
 - Attach the mixer tip, aligning the projections on the mixer tip with the grooves on the syringe, and turn 90° in a clockwise direction. (Fig. C, Fig. D)
 - For post cementation
 - Set the mixer tip and nozzle for endodontic treatment following the standard method.
4. Replace the cap securely after use.
 - Replace the mixer tip for each use.
 - Thoroughly remove the paste remaining inside the cap before replacing it.

Preparation of handmix syringe

1. Uncap the handmix syringe.
 - Do not remove the syringe nozzle.
2. Gently place the nozzle tips onto a mixing pad and depress the plunger to dispense the pastes, confirming that the pastes are flowing evenly at the same level.
3. After dispensing, clean the nozzle tips with a gauze etc., and replace the cap securely.
 - Make sure the right direction of the cap before replacing it. (Fig. E)


Fig. E

A. Intraoral cementation of restorations

A-1. Preparation of the restoration

Roughen the adhesive surface of the restoration by etching, using abrasives, or sandblasting with 50-100 µm alumina particles according to Table 1. Rinse with water and dry.

(Table 1) Adhesive surface	Pressure
Metal	Approx. 0.3-0.5 MPa (Approx. 3-5 kgf/cm ²)
Alumina, zirconia	Approx. 0.2-0.3 MPa (Approx. 2-3 kgf/cm ²)
Silicate glass, composite resin	Approx. 0.1-0.2 MPa (Approx. 1-2 kgf/cm ²)
CAD/CAM hybrid ceramics	Approx. 0.2-0.3 MPa (Approx. 2-3 kgf/cm ²)

- Perform appropriate pretreatment and apply optimum pressure following the manufacturers' instructions for use of each material.

A-2. Preparation of cavity or abutment

For natural tooth abutment

Remove all the temporary materials thoroughly. Clean the cavity or abutment tooth surface and isolate the treatment area following the standard method. Protect the pulp if necessary.

- In the area where close proximity to the pulp, apply calcium hydroxide to protect the pulp or use glass ionomer/resin-modified glass ionomer cement as a base/liner.

For metal or resin abutment

Prepare and clean the abutment following the standard method.

✓ *To enhance adhesion, pretreat the adhesive surfaces of the abutment/ cavity and prosthesis with BeautiBond Xtreme following their instructions for use.*

A-3. Application of resin cement

Automix syringe

Gently press the plunger to apply BeautiLink SA to the adhesive surface of the restoration. Handmix syringe
Depress the plunger to dispense the pastes and mix them sufficiently using a plastic spatula for 10 sec. Apply BeautiLink SA to the adhesive surface of the restoration.

- Do not apply BeautiLink SA directly to the tooth surface in patient's mouth.

A-4. Seating of the restoration

Seat the restoration in the cavity or onto the abutment and apply proper pressure.

- Using polyester film etc., between preparation and adjacent teeth prior to placement of the restoration aids in isolation and easy removal of excess cement.
- A-5. Removal of excess cement**
When using a light-curing unit
Light-cure the excess material for 1-2 sec with a light-curing unit (Halogen or LED). Then, remove the tack-cured material with a probe or other suitable instrument.
For the area where light does not reach; or when using Opaque cement
Leave 2-3 min after placing the restoration to allow the material to cure chemically. Then, remove the excess material with a probe or other suitable instrument. Or remove it with a cotton ball before self-curing completely.
➢ When using protective material (e.g., matrix band), remove it within 2 min from start of mixing the paste.

A-6. Curing of cement

Light-cure each surface and marginal area with a light-curing unit; Halogen for 20 sec/LED for 10 sec. However, for the area where the light does not reach, or when using Opaque cement, hold the restoration in place for around 5 min after light-curing in order to allow it to cure chemically.
➢ Light-cure the margin sufficiently. Sufficient light-curing of marginal area provides more stable cementation of CAD/CAM hybrid ceramics etc.

A-7. Finishing

If necessary, polish the margin and adjust occlusion following the standard method.

B. Extraoral cementation of the restorations

B-1. Preparation of the restoration

Roughen the adhesive surface of the restoration by etching, using abrasives, or sandblasting with 50-100 µm alumina particles according to Table 1. Rinse with water and dry.

B-2. Application of resin cement

Automix syringe

Apply BeautiLink SA to the adhesive surface of the restoration. Handmix syringe
Depress the plunger to dispense the pastes and mix them sufficiently using a plastic spatula for 10 sec. Apply BeautiLink SA to the adhesive surface of the restoration.

B-3. Cementation of the restorations

Cement the restorations and hold them together with proper pressure.

B-4. Removal of excess cement

Leave 3-4 min after cementing the restorations to allow the material to cure chemically. Then, remove excess material with a probe or other suitable instrument. Or remove it with a cotton ball before self-curing completely.

B-5. Curing of the cement

Hand held light-curing unit
Light-cure each surface and marginal area with a light-curing unit; Halogen for 20 sec/LED for 10 sec.
Lab light-curing unit
Light-cure the material for 180 sec.
➢ When cementing a restoration that does not transmit light, such as metal, or when using Opaque cement, allow the material to self-cure by holding the restoration for approx. 10 min after light-curing.
➢ Do not heat-cure the material.

B-6. Finishing

If necessary, polish the margin following the standard method.

C. Post cementation

- ✓ *Opaque shade is not recommended for post cementation.*

C-1. Preparation of post space

Prepare and fill the root canal following the standard method. Then, prepare the post space. Isolate with a rubber dam etc., if necessary.

C-2. Preparation of the post (trial seat)

Choose an adequate size (length and diameter) of post and try seating it into the post space. Clean the post with alcohol moistened cotton ball etc., and treat the surface of the post properly following manufacturer's instructions for use.

C-3. Pretreatment of the post space

Apply BeautiBond Xtreme sufficiently onto the entire adhesive surface. Air dry with gentle air for 3 sec and then, dry with stronger air to dry the surface completely. Light-cure with a light-curing unit; Halogen for 10 sec/LED for 5 sec.
➢ If the adhesive pools in the root canal, absorb excess material with a paper point etc. Then air dry.
➢ Irradiate as closely as possible to the restoration area. In cases with a deep post space, additional light-curing is recommended.

C-4. Post cementation

Fill BeautiLink SA into the post space. Insert the post into the post space immediately after filling the material. Then, light-cure from every aspect using a light-curing unit to fix the post; Halogen for 20 sec/LED for 10 sec.

C-5. Core build-up

Build up the core using a core material following the standard method and prepare the abutment.

(Notes)

- Use hand held light-curing unit with the specifications below.

Light-curing unit	Halogen	LED
Wavelength	400-500 nm	440-490 nm
Light intensity	≥ 500 mW/cm ²	≥ 1,000 mW/cm ²

- Operators must use protective glasses or light shielding plate to avoid direct exposure to the curing light. Operators must also protect patients' eyes from the direct exposure to the curing light.
- Disinfect the light-curing unit with alcohol after each use to avoid cross contamination.
- Do not use eugenol-containing material and hydrogen peroxide since they might retard the curing process and weaken the bond strength.
- Use saliva ejector or vacuum during treatment.
- Avoid contamination of the adhesive surface by saliva or blood, etc., before bonding. If the surface is contaminated by saliva or blood, clean thoroughly with water or an ethanol moistened cotton ball and dry. Then, reapply the material.
- Remove any hardened paste plugging the opening of the syringe before use.
- Mixer tips are for single use only.
- Thoroughly remove the paste remaining inside the cap before replacing it. Make sure the cap is securely closed to avoid polymerization due to ambient light.
- Activation by external energy is mandatory.
- Follow the instructions for use of each dental material, instrument, or equipment to be used in conjunction with this product.
- Do not use this product for any purposes other than specifically outlined in these instructions for use.
- This product is intended for use by dental professionals only.
- After opening the foil pouch, use this product within 6 months and before the expiration date.
- When you open the pouch, make sure to fill out the opening date on the sticker on the pouch to put it on the syringe.

POSSIBLE SIDE EFFECTS/RESIDUAL RISKS

- This product contains substances that may cause allergic reactions.
- Debonding and poor marginal adaptation of the restoration due to long-term use might require re-treatment.

COMPOSITION

Paste A	
Composition	% by weight
Zirconium silicate	40-60
Bis-GMA	10-20

Composition	
	% by weight
Dimethacrylate monomer	10-20
Acid monomer	5-20
Silicate glass	1-10
Polymerization initiator	< 5
Others	< 5

Paste B

Composition	
	% by weight
Fluoroboroaluminosilicate glass	40-60
UDMA	20-30
Silane coupling agent	5-15
Dimethacrylate monomer	1-10
Silicate glass	1-10
Polymerization initiator	< 5
Others	< 5

SHADES

Clear, Ivory, Opaque

STORAGE

- Store at 1-25 °C/34-77 °F.
- Avoid high temperatures and high humidity.
- Keep away from direct sunlight and any source of ignition.

DISPOSAL

Dispose the contents and containers in accordance with relevant local and national regulations.

PHYSICAL PROPERTIES

Setting time (37 °C/98 °F)	≤ 8 min
Working time (23 °C/73 °F)	≥ 110 sec
Radio-opacity	≥ 1.25 mm of Al*
Bond strength to dentin	≥ 5 MPa

In accordance with ISO/TS 16506

*Presents 1.25 times or higher radio-opacity than natural dentin.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

Location of the symbols glossary: https://www.shofu.com/en/symbols_glossary/

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