

Beautiful-Bulk Flowable

Bulk Fill Flowable Composite with Bioactive GIOMER Technology




Beautiful-Bulk Flowable

Bulk Fill Flowable Composite with Bioactive GIOMER Technology


Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com
800.827.4638
customer-service@shofu.com




FUNCTIONAL

Beautiful-Bulk Flowable is the ideal light-cure dentin replacement material for deep restorations. For years glass ionomer was used for dentin replacement, but long curing times caused significantly reduced productivity, making this impractical until resin-based alternatives emerged. With physical properties similar to dentin, 4mm depth of cure, and the lowest shrinkage stress in this category, **Beautiful-Bulk Flowable** is the perfect choice for sensitivity-free dentin replacement and maximum clinical efficiency.




BIOACTIVE GIOMER TECHNOLOGY INSIDE

Shofu's proprietary Giomer Technology gives **Beautiful-Bulk Flowable** anti-bacterial and acid neutralizing effects to minimize problems caused by bacterial proliferation such as secondary decay and restorative failure. A Univ. of Florida Dental School clinical study showed no secondary decay or failures after 8 years! **Beautiful-Bulk Flowable** is the only bulk flow product to combine bioactive benefits with speed and clinical efficiency.




Minimal Bacteria (Asahi University)




Advanced Bacterial Colonization

Fluoride Release & Recharge from GIOMER



*1. Gordan V.V., et al. A clinical evaluation of a self-etching primer and a Giomer restorative material: Results at 8 years. JAM Dent. Assoc. 2007;138:621-627

Beautiful-Bulk Flowable



Beautiful-Bulk Flowable

2.4 g syringe

Shade

Universal

Dentin

Part Number

2030

2031

0.23 g x 20 tips

Shade

Universal

Dentin

Part Number

2028

2029



Scan Here
Instructions for Use