

# Beautiful-Bulk Restorative

## Bulk Fill Posterior Composite with Bioactive GIOMER Technology




### Beautiful-Bulk Restorative

Bulk Fill Posterior Composite with Bioactive GIOMER Technology

Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com  
800.827.4638  
customer-service@shofu.com



### FUNCTIONAL


**Beautiful-Bulk Restorative** is the ideal metal-free material for bulk fill posterior restorations. This highly filled material (87.0%) offers low shrinkage (<1.7%) and wear resistance, 4mm depth of cure, as well as the combination of superior compressive and flexural strength for maximum durability. In addition, **Beautiful-Bulk Restorative's** esthetics are unmatched by any competitor.




Jack Griffin, DDS

### BIOACTIVE GIOMER TECHNOLOGY INSIDE

Shofu's proprietary GIOMER Technology gives **Beautiful-Bulk Restorative** long-term ion exchange for anti-bacterial and acid neutralizing effects to minimize problems caused by plaque and biofilm which may lead to restorative failure. A Univ. of Florida Dental School clinical study showed no secondary decay or failures after 8 years! **Beautiful-Bulk Restorative** is the only product in this market segment that gives the bioactive and therapeutic benefits of GIOMER Technology and the clinically efficiency of posterior bulk filling.




Minimal Bacteria (Asahi University)




Advanced Bacterial Colonization

#### Fluoride Release & Recharge from GIOMER



\*1. Gordan V.V., et al. A clinical evaluation of a self-etching primer and a Giomer restorative material: Results at 8 years. JAM Dent. Assoc. 2007;138:621-627

## Beautiful-Bulk Restorative



### Beautiful-Bulk Restorative

4.5 g syringe

Shade  
Universal  
A

Part Number

2034  
 2035

0.25 g x 20 tips

Shade  
Universal  
A

Part Number

2032  
 2033



Scan Here

Instructions for Use