














SHOFU ROBOT® CARBIDE CTF-FG FINISHERS



- High quality tungsten carbide finishing head on a stainless steel shank
- Trim & finish subgingival cosmetic restorations without damaging the subgingival tissue
- Available in three blade types: 8, 16 (yellow-banded shank) & 30 (white-banded shank)
- 8 blades: rapid trimming & contouring; 16 blades: Fine Finishing (F); 30 blades: Ultra Fine Finishing (UF)
- **Use Robot Carbide CTF-FG Finishers for:** cosmetic restorations

SHAPE CATEGORY		SH 132 F	SH 132 UF	SH 133	SH 133F	SH 133 UF
SHAPE						
HEAD SIZE (MM)		008	008	010	010	010
HEAD LENGTH (MM)		3.0	3.0	4.2	4.2	4.2
# OF BLADES		16	30	8	16	30
MAXIMUM RPM		300,000	300,000	450,000	450,000	450,000
FG (3 PK)	PN	0762-3	0763-3	0764-3	0765-3	0766-3

SHAPE CATEGORY		SH 134	SH 134 F	SH 134 UF	SH 135	SH 135F	SH 135 UF
SHAPE							
HEAD SIZE (MM)		014	014	014	014	014	014
HEAD LENGTH (MM)		6.0	6.0	6.0	9.0	9.0	9.0
# OF BLADES		8	16	30	8	16	30
MAXIMUM RPM		450,000	450,000	450,000	300,000	300,000	300,000
FG (3 PK)	PN	0767-3	0768-3	0769-3	0770-3	0771-3	0772-3

PRICING - SHOFU ROBOT® CARBIDE CTF-FG FINISHERS

<p>PN</p> <p>See charts CTF-FG Regular - 8 blades (3 per pack)</p> <p>See charts CTF-FG Fine - 16 blades (3 per pack)</p> <p>See charts CTF-FG Ultrafine - 30 blades (3 per pack)</p>	<p>Product</p>
--	-----------------------