

WINDOWS

Program	EYECAM
Adstra	Y
Apixia Digirex	U
Ateco	U
Choice/CADI	Y
Cleardent	Y**
ClioSoft	N
Curve Hero	Y#
Dental Eye/MiPACS	Y
DentalEMR	U
Dentrix Image	Y
Dentrix Ascend	Y**
Dexis	Y**
Eaglesoft/Patterson Imaging	Y**
Easydental	N
Easydent4	Y
Eva Prolmage	U
EvaSoft	U
EZDent-i	Y
EZImage	U
Genesis	Y
IdentalSoft	YT
Kodak Dental Imaging	Y**
Mediadent V5	U
Mogo	N
MyRay Raymage	N
Opendental	Y
Oryx	Y
Owandy Quick Vision	N
Patient Gallery	Y
Progeny Vet Imaging	U
QuickImage DIS (Archimed)	Y
Romexis	Y
Safecom	Y
Scanora	Y**
Schick CDR Dicom	N
Sidexis	N
SoftDent Imaging	N
Sopro Imaging	Y**

WINDOWS (CONTINUED)

Program	EYECAM
Sota Imaging	YT
Suni	Y**
Tigerview/Quickview	Y
Tracker Imaging	U
Vatech EasyDent 4	U
Vatech EasyDent-i	U
Video Dental Imaging	Y**
Visix \$	N
VixWin	Y
XDR	Y**
XV/XV Lite	Y**

APPLE MAC OS

Program	EYECAM
MacPractice	Y
Curve Hero	Y
Dentrix Ascend	Y
RadioVision	Y
Dovetail(web)	U

Legend

Y#	Yes, but requires Curve to be run in Chrome browser
U	Untested
Y	Yes
Y*	Yes, but requires DentCapture software to capture
Y**	Yes, but requires DentCapture software to capture
N	No
N?	No But Needs To Be Verified
\$	Requires fee to be paid to software manufacturer to work
YT	Yes, works via TWAIN
T.O.	Only can import from another program with twain to capture
YF	Yes, but requires foot pedal to capture
YV	Yes, but requires the use of VideoGlide