

OneGloss PS

One System - One Step
Finishing and Polishing



OneGloss PS
One System - One Step
Finishing and Polishing

Shofu is dedicated to manufacturing the highest quality products utilizing our innovative technologies to enhance minimally invasive cosmetic dentistry.

shofu.com
800.827.4638
customer-service@shofu.com



FAST AND EASY

Save chairside time when finishing and polishing all composite restorations with **OneGloss PS**. When used on the slow speed with water spray; a firmer pressure will remove composite while softer pressure will polish to a satin finish.



Finish and polish in one step.



Lasts longer than the competitors.




Does not scratch enamel.

PROVEN QUALITY AND SUPERIOR VALUE



When compared to Enhance®, the leading competitor, **OneGloss PS** is rated much higher; finishing and polishing more restorations before showing any signs of deterioration.* Even better, **OneGloss PS** offers 25% more product for about 20% less.†



†Information is based on Henry Schein Dental Catalog.

OneGloss

PROVEN QUALITY

For almost 100 years, Shofu Dental's dedication to quality and value have made us one of the world's most respected dental products manufacturers.

SHAPE CATEGORY	Cup	Midi-Point	IC	Assort. Kit
Shape				-
OneGloss				
CA (50 pk) PN	0181	0182	0183	0180
OneGloss PS				
CA (50 pk) PN	0176	0177	0178	0175